

Choose the frame and seat combination to find the max. user load

R82[®]

| Product | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|----------------------------|-----------------------------|----------------------------|----------------------------|----------------------------|------------------------------|------------------------------|-------------------------------|------------------------------|------------------------------|------------------------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|--------------------|--------------------|-------------------|-------------------|-------------------|-------------------|----------------------|----------------------|---|---|---|---|---|
| | Panda Futura Active Size 2 | Panda Futura Active Size 2½ | Panda Futura Active Size 3 | Panda Futura Active Size 4 | Panda Futura Active Size 5 | Panda Futura Standard Size 1 | Panda Futura Standard Size 2 | Panda Futura Standard Size 2½ | Panda Futura Standard Size 3 | Panda Futura Standard Size 4 | Panda Futura Standard Size 5 | x:panda Size 1S | x:panda Size 1M | x:panda Size 1L | x:panda Size 2S | x:panda Size 2M | x:panda Size 2L | x:panda Size 3S | x:panda Size 3M | x:panda Size 3L | x:panda Size 4S | x:panda Size 4M | x:panda Size 4L | Serval seat Size 1 | Serval seat Size 2 | Nandu seat Size 1 | Nandu seat Size 2 | Nandu seat Size 3 | Nandu seat Size 4 | Stingray seat Size 1 | Stingray seat Size 2 | | | | | |
| High-Low Indoor Size 1 | Ordinary use | 50 | 50 | 50 | 50 | - | 50 | 50 | 50 | 50 | 50 | 45* | 45* | 45* | 50* | 50* | 50* | - | - | - | - | - | 35 | 43 | 50 | 50 | 50 | - | - | | | | | | | |
| | Ordinary use | 50 | 50 | 50 | 50 | - | 50 | 50 | 50 | 50 | 50 | 45* | 45* | 45* | 50* | 50* | 50* | - | - | - | - | - | 35 | 43 | - | - | - | - | - | - | | | | | | |
| High-Low Outdoor Size 1 | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 27 | 43 | 43 | - | 43* | 43* | 43* | 43* | 43* | 43* | - | - | - | - | - | - | - | - | - | - | - | - | | | | | |
| High-Low:x Indoor Size 1 | Ordinary use | 50* | 50* | 50* | 50* | - | 50* | 50* | 50* | 50* | 50* | 45 | 45 | 45 | 50 | 50 | 50 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | | | | | |
| High-Low:x Indoor Size 2 | Ordinary use | 50* | 50* | 50* | 50* | 43 | 50* | 50* | 50* | 50* | 50* | 43 | - | - | - | - | - | - | 47 | 47 | 47 | - | - | - | - | - | - | - | - | - | - | | | | | |
| High-Low:x Indoor Size 3 | Ordinary use | - | - | - | 50* | 85 | - | - | - | 50* | 83 | - | - | - | - | - | - | 60 | 60 | 60 | 80 | 80 | 80 | - | - | - | - | - | - | - | - | - | - | | | |
| | Ordinary use | 50* | 50* | 50* | 50* | - | 50* | 50* | 50* | 50* | 50* | 45 | 45 | 45 | 50 | 50 | 50 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | | | | |
| High-Low:x Outdoor Size 1 | Transportation in vehicles | - | - | - | - | - | - | 27* | 27* | 27* | 43* | 43* | - | 43 | 43 | 43 | 43 | 43 | 43 | - | - | - | - | - | - | - | - | - | - | - | - | - | | | | |
| | Ordinary use | 50* | 50* | 50* | 50* | 43 | 50* | 50* | 50* | 50* | 50* | 43 | - | - | - | - | - | - | 60 | 60 | 60 | - | - | - | - | - | - | - | - | - | - | - | - | | | |
| High-Low:x Outdoor Size 2 | Transportation in vehicles | - | - | - | - | - | 43 | 27* | 27* | 27* | 43* | 43* | - | 43 | - | - | - | - | 43 | 43 | 43 | - | - | - | - | - | - | - | - | - | - | - | - | | | |
| High-Low:v Indoor Size 1 | Ordinary use | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | 40 | - | - | | | |
| High-Low:v Indoor Size 2 | Ordinary use | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | 40 | 40 | - | - | | | |
| Nandu Frame Size 1 | Ordinary use | 50* | 50* | - | - | - | 50* | 50* | 50* | - | - | 45* | 45* | 45* | 50* | 50* | 50* | - | - | - | - | - | - | - | - | 50 | 50 | - | - | - | - | - | - | - | - | |
| Nandu Frame Size 2 | Ordinary use | - | - | 50* | 50* | - | - | - | 50* | 50* | - | - | - | 50* | 50* | 50* | 50* | 50* | 50* | - | - | - | - | - | - | - | 50 | 50 | - | - | - | - | - | - | - | - |
| Combi Frame:x Size 1 20" 12½" | Ordinary use | 50 | 50 | 50 | - | - | 50 | 50 | 50 | 50 | 50 | 45 | 45 | 45 | 50 | 50 | 50 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 43 | 43 | - | 45 | 45 | 45 | 50 | 50 | 50 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | |
| Combi Frame:x Size 1 22" 24" | Ordinary use | 50 | 50 | 50 | - | - | 50 | 50 | 50 | 50 | 50 | 45 | 45 | - | 50 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 43 | 43 | - | 45 | 45 | - | 50 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | |
| Combi Frame:x Size 2 20" 22" 24" 12½" | Ordinary use | 50 | 50 | 50 | 50 | 58 | 50 | 50 | 50 | 50 | 50 | 45 | 45 | 45 | 50 | 50 | 50 | 60 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 43 | 43 | - | 45 | 45 | 45 | 50 | 50 | 50 | 57 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | |
| Combi Frame:x Size 3 20" 12½" | Ordinary use | 50 | 50 | 50 | 50 | 58 | 50 | 50 | 50 | 50 | 50 | 45 | 45 | 45 | 50 | 50 | 50 | 60 | 60 | 60 | 60 | - | - | - | - | - | - | - | - | - | - | - | - | | | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 43 | 43 | - | 45 | 45 | 45 | 50 | 50 | 50 | 57 | 57 | 57 | - | - | - | - | - | - | - | - | - | - | - | - | - | | |
| Combi Frame:x Size 3 22" 24" | Ordinary use | 50 | 50 | 50 | 50 | 58 | 50 | 50 | 50 | 50 | 50 | 45 | 45 | 45 | 50 | 50 | 50 | 60 | 60 | 60 | 60 | - | - | - | - | - | - | - | - | - | - | - | - | | | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 43 | 43 | - | 45 | 45 | 45 | 50 | 50 | 50 | 57 | 57 | 57 | - | - | - | - | - | - | - | - | - | - | - | - | - | | |
| Combi Frame:x Size 4 20" 22" 24" 12½" | Ordinary use | 50 | 50 | 50 | 50 | 83 | 50 | 50 | 50 | 50 | 50 | 45 | 45 | 45 | 50 | 50 | 50 | 60 | 60 | 60 | 80 | 80 | - | - | - | - | - | - | - | - | - | - | - | | | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 43 | 43 | 85 | 45 | 45 | 45 | 50 | 50 | 50 | 57 | 57 | 57 | 57 | - | - | - | - | - | - | - | - | - | - | - | - | | |
| Multi Frame Size 1 - 12½" | Ordinary use | 50 | 50 | 50 | - | - | 50 | 50 | 50 | 50 | 50 | 45 | 45 | 45 | 50 | 50 | 50 | 60 | 60 | 60 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 27 | 50 | - | 45 | 45 | 45 | 50 | 50 | 50 | 57 | 57 | 57 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | |
| Multi Frame Size 1 - 22" | Ordinary use | 50 | 50 | 50 | - | - | 50 | 50 | 50 | 50 | 50 | 45 | 45 | - | 50 | 50 | 50 | 60 | 60 | 60 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | |
| | Transportation in vehicles | - | - | - | - | - | - | 27 | 27 | 27 | 50 | - | 45 | 45 | - | 50 | 50 | 50 | 57 | 57 | 57 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | |
| Multi Frame Size 2 - 12½" | Ord | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |