

# Important markings

Before cleaning your x:panda shape, there are some important markings that can be helpful to note down. In this way you can easily get back to the needed settings.



## Wedges



	1	1½	2	2½	3	3½	4	4½	5	5½	6	6½	7	7½	8	8½	9	9½	10	10½	11	11½	12	12½	13
Users right																									
Users left																									



## Inner and outer knee supports

	S	M	L		S
Inner knee supports					Outer knee supports



## Seat plates for advanced seating

Users right							Users left						
3	2	1	0	1	2	3	3	2	1	0	1	2	3




	10	20	30	40	50	60	70	80	90	100	110	120	130	140
Users right														
Users left														





Explore our how-to-videos to learn more about the adjustments on the x:panda shape seating system.

# Important markings

## Fix locks at the lower back and upper back for mounting of vest



## Fix locks at the seat base for mounting of belt



Explore our how-to-videos to learn more about the adjustments on the x:panda shape seating system.