

# Choose the frame and seat combination to find the max. load



Product		Panda Futura Active Size 2	Panda Futura Active Size 2½	Panda Futura Active Size 3	Panda Futura Active Size 4	Panda Futura Active Size 5	Panda Futura Standard Size 1	Panda Futura Standard Size 2	Panda Futura Standard Size 2½	Panda Futura Standard Size 3	Panda Futura Standard Size 4	Panda Futura Standard Size 5	x.panda Size 1S	x.panda Size 1M	x.panda Size 1L	x.panda Size 2S	x.panda Size 2M	x.panda Size 2L	x.panda Size 3S	x.panda Size 3M	x.panda Size 3L	x.panda Size 4S	x.panda Size 4M	x.panda Size 4L	Stingray seat Size 1	Stingray seat Size 2
High-Low Indoor Size 1	Ordinary use	30	50	50	50	-	30	30	50	50	50	-	30*	30*	30*	30*	30*	30*	-	-	-	-	-	-	-	-
High-Low Outdoor Size 1	Ordinary use	30	50	50	50	-	30	30	50	50	50	-	30*	30*	30*	30*	30*	30*	-	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	-	30	30	43	43	43	-	30*	30*	30*	30*	30*	30*	-	-	-	-	-	-	-	-
High-Low:x Indoor Size 1	Ordinary use	30	50*	50*	50*	-	30*	30*	50*	50*	50*	-	30*	30*	30*	30*	30*	30*	-	-	-	-	-	-	-	-
High-Low:x Indoor Size 2	Ordinary use	30	50*	50*	50*	43	30*	30*	50*	50*	50*	60	-	-	-	-	-	-	50	50	50	-	-	-	-	-
High-Low:x Indoor Size 3	Ordinary use	-	-	-	60	85	-	-	-	-	60*	85	-	-	-	-	-	-	60	60	60	80	80	80	-	-
High-Low:x Outdoor Size 1	Ordinary use	30	50*	50*	50*	-	30*	30*	50*	50*	50*	-	30	30	30	30	30	30	-	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	-	30*	30*	43*	43*	43*	-	30	30	30	30	30	30	-	-	-	-	-	-	-	-
High-Low:x Outdoor Size 2	Ordinary use	30	50*	50*	50*	43	30*	30*	50*	50*	50*	60	-	-	-	-	-	-	50	50	50	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	43	30*	30*	43*	43*	43*	43	-	-	-	-	-	-	43	43	43	-	-	-	-	-
High-Low:v Indoor Size 1	Ordinary use	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	40	-
High-Low:v Indoor Size 2	Ordinary use	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	40	40
Combi Frame:x Size 1 12½"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	60	80	-	-	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	60	75	-	-	-	-
Combi Frame:x Size 1 20" 22" 24"	Ordinary use	30	50	50	-	-	30	30	50	50	-	-	30	30	30	30	30	-	-	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	-	30	30	50	50	-	-	30	30	30	30	30	-	-	-	-	-	-	-	-	-
Combi Frame:x Size 2 12½"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	60	80	80	-	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	60	75	75	-	-	-
Combi Frame:x Size 2 20" 22" 24"	Ordinary use	30	50	50	60		30	30	50	50	60		30	30	30	30	30	30	60	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-		30	30	50	50	57		30	30	30	30	30	30	60	-	-	-	-	-	-	-
Combi Frame:x Size 3 12½"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	60	80	80	80	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	60	75	75	75	-	-
Combi Frame:x Size 3 20" 22" 24"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	-	-	-	-	-	-
Combi Frame:x Size 4 12½"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	60	80	80	80	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	60	75	75	75	-	-
Combi Frame:x Size 4 20" 22" 24"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	60	80	80	-	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	60	75	75	-	-	-
Multi Frame Size 1 - 12½"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	-	-	-	-	-	-
Multi Frame Size 1 - 22"	Ordinary use	30	50	50	-	-	30	30	50	50	-	-	30	30	-	30			-	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	-	30	30	50	50	-	-	30	30	-	30			-	-	-	-	-	-	-	-
Multi Frame Size 2 - 12½"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	60	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	60	-	-	-	-	-
Multi Frame Size 2 - 22"	Ordinary use	30	50	50	60	-	30	30	50	50	60	-	30	30	30	30	30	30	60	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	-	30	30	50	50	57	-	30	30	30	30	30	30	60	-	-	-	-	-	-	-
Multi Frame Size 3 - 12½"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	60	80	80	80	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	60	75	75	75	-	-
Multi Frame Size 3 - 22"	Ordinary use	30	50	50	60	85	30	30	50	50	60	85	30	30	30	30	30	30	60	60	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	85	30	30	50	50	57	85	30	30	30	30	30	30	60	60	-	-	-	-	-	-
Multi Frame:x Size 3 - 12½"	Ordinary use	-	-	-	-	85	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	80	80	80	-	-
	Transportation in vehicles	-	-	-	-	85	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	75	75	75	-	-
Multi Frame:x Size 3 - 22"	Ordinary use	-	-	-	-	85	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	85	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Stingray Frame	Ordinary use	-	-	-	-	-	30	30	40	-	-	-	30	30	30	30	30	-	-	-	-	-	-	-	40	40
	Transportation in vehicles	-	-	-	-	-	30	30	40	-	-	-	30	30	30	30	30	-	-	-	-	-	-	-	40	40
Strong Base	Ordinary use	-	-	-	-	85	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	80	80	80	-	-

\* Special interface needed