

Ropox Product Improvements

Two of our changing beds just got even better: The **ROPOX Shower/Changing Bed** and the **ROPOX Vario Changing Bed**. And it's not just a tiny little thing we've adjusted. We've made the changing beds better for everyone: the user, care worker and service engineer.



Ropox Shower/Changing Bed



Ropox Vario Changing Bed



New hinge mechanism



Old hinge mechanism

Product improvements

The following construction improvements apply to both changing bed models:

- New redesigned hinge mechanism eliminating sharp edges – pls. see pictures above
- More slim, simple and discreet unit (depth before: 127mm. / now: 75mm.) – pls. see picture above
- Easier to remove cover making service of electrical components simpler
- Better possibility of fixing the cover hereby reducing the risk of vandalism
- Power outlet can now be placed behind the unit
- Easier installation and service in low-ceiling rooms (under 234cm.) - The cover can now be removed without pulling the cover up and over the top of the unit

Some things haven't changed

The model types, accessories and item numbers are the same.

The products' key elements also remain unchanged: the long electrical height adjustment range, the reversible, soft slats, and the ability to fold up the changing bed and save space. All the things that make them safe and comfortable for the user and ergonomic and practical for the care workers.